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ENGLISH

DESCRIPTION

GalaFLEX 3DR™ scaffold is a single layer, bioresorbable surgical mesh made from poly-4-hydroxybutyrate (P4HB). P4HB is produced from a naturally occurring monomer and is processed into monofilament fibres and knitted into a surgical scaffold. The scaffold has a 3D curvature with a supporting rim designed to promote better conformance with a patient's anatomy in locations in which a flat design does not easily conform. P4HB bioresorbs through a process of hydrolysis and hydrolytic enzymatic digestion. It has been developed to optimise resorption rate and prolong strength retention in order to provide support throughout the expected period of healing. Although the scaffold loses strength with time, its porous construction was designed to allow native tissue ingrowth and gradual transfer of load from the scaffold to the tissue.

INTENDED USE

GalaFLEX 3DR™ scaffold is intended for use, as an adjunct to sutures, for the reinforcement and repair of soft tissue where weakness exists and where the addition of a reinforcing material is needed to obtain the desired surgical result in patients undergoing breast surgery. The GalaFLEX 3DR™ scaffold is designed to be used in patients undergoing soft tissue repair and reinforcement in medically necessary breast surgery procedures where the existing soft tissue is deficient to support the surgical repair. Examples of such breast surgery applications include reduction mammoplasty and breast revision surgery to correct a medical condition. GalaFLEX 3DR™ scaffold may also be used in cosmetic breast procedures.

CONTRAINDICATIONS

None known.

WARNINGS

1. Device manufacture involves exposure to tetracycline hydrochloride and kanamycin sulfate. The safety and product use for patients with hypersensitivities to these antibiotics is unknown.
2. The safety and effectiveness of GalaFLEX 3DR™ scaffold in neural tissue and in cardiovascular tissue has not been established.
3. The safety and effectiveness of GalaFLEX 3DR™ scaffold in paediatric use has not been established.
4. Placement of the scaffold in direct contact with bowel or viscera is not recommended.
5. If an infection develops, treat the infection aggressively. An unresolved infection may require removal of the scaffold.
6. Because GalaFLEX 3DR™ scaffold is fully bioresorbable, it should not be used in repairs where permanent support from the scaffold is required.
7. GalaFLEX 3DR™ scaffold is supplied sterile. Inspect the device and packaging prior to use to be sure they are intact and undamaged.
8. GalaFLEX 3DR™ is for single use only. Do not re-sterilise or re-use any portion of GalaFLEX 3DR™ scaffold. Never re-use a scaffold, even if it

Printed on: 10 Apr 2024, 11:37:49 am; Printed by: KEARLE.

seems to be in perfect condition, in order to prevent risk of cross-contamination or risk of reduced performance. Do not re-sterilise the scaffold. Once the package has been opened, any unused scaffold should be discarded.

9. Unused scaffold must be discarded according to the institution's procedures for handling of biohazardous materials.
10. Any decision to remove the scaffold should take into account potential risks associated with a second surgical procedure. This may include difficulty to remove the scaffold due to the ingrowth of tissue or degradation of the scaffold. Scaffold removal should be followed by adequate post-operative management.
11. The GalaFLEX 3DR™ scaffold has not been studied for use in breast reconstructive surgeries.

PRECAUTIONS

Only doctors qualified in the appropriate surgical techniques should use this device. Users should be familiar with surgical procedures and techniques, including strength requirements and scaffold size choices. Improper selection, placement, positioning and fixation of GalaFLEX 3DR™ can cause subsequent undesirable results. Patient should be provided with instructions regarding post-operative care, i.e. lifting, hygiene, activity limitations and any other specific patient/procedure requirements. GalaFLEX 3DR™ scaffold is bioresorbable and therefore any initial palpability of the scaffold decreases over time.

The safety and effectiveness of GalaFLEX 3DR™ scaffold in the proximity of existing or excised cancer has not been established.

ACTIONS

GalaFLEX 3DR™ scaffold degrades through a process of hydrolysis and hydrolytic enzymatic digestion.

It has been developed to minimise the variability of resorption rate and strength and provide support throughout the expected period of healing. Pre-clinical implantation studies indicate that GalaFLEX 3DR™ scaffold retains approximately 70% of its strength at 12 weeks. Absorption of the scaffold material will be essentially complete within 18-24 months. The outer rim may be palpable throughout the expected period of healing.

ADVERSE REACTIONS

Possible complications of using GalaFLEX 3DR™ scaffold include infection, seroma, pain or swelling, scaffold migration, wound dehiscence, haemorrhage, adhesions, haematoma, inflammation, extrusion and recurrence of the soft tissue defect. In pre-clinical testing, GalaFLEX 3DR™ scaffold elicited a minimal tissue reaction characteristic of foreign body response to a substance. The tissue reaction resolved as the scaffold was resorbed.

IMAGING COMPATIBILITY

GalaFLEX 3DR™ scaffold is made from fully resorbable P4HB and is temporary in nature. This material is non-conducting, non-metallic and non-magnetic. Therefore, in accordance with the definition stated in ASTM F-2503-13, Standard Practice for Marking Devices and Other Items for Safety in the Magnetic Resonance
Printed on: 10 Apr 2024, 11:37:49 am; Printed by: KEARLE.

Environment, the device is determined to be "MR Safe - an item that poses no known hazards in all MR environments."

DIRECTIONS FOR USE

1. Prior to implanting GalaFLEX 3DR™ scaffold, complete repair using appropriate suture technique.
2. GalaFLEX 3DR™ scaffold should be prepared for placement using standard surgical preparation technique.
3. GalaFLEX 3DR™ scaffold is not intended to be cut as the rim enables the 3D curvature to be maintained. If the rim is cut, the scaffold may lose its three dimensional shape and may not conform as readily to the anatomical site. Using aseptic technique, fixation tabs may be cut if desired for each specific application.
4. To ensure the scaffold remains curved, minimise pulling and stretching of the scaffold prior to implantation.
5. Transfer the scaffold to the surgical site. Implant the scaffold so that the contour of the scaffold follows the natural contours of the implantation site and the edge of the scaffold extends beyond the margins of the defect. Suture the scaffold into place, avoiding excessive tension and without the expectation of stretch.
6. It is recommended that suture fixation be placed 6 mm to 12 mm apart at a distance approximately 6 mm from the edge of the scaffold using an interrupted suturing technique.
7. If using a resorbable suture to fixate the scaffold, please ensure the strength retention of the suture is appropriate for your application. The edges or corners of the scaffold should be fixated with suture such that it lies flat against the tissue of the repair site. The scaffold should be sufficiently anchored to assure proper closure under correct tension and to stabilise it during tissue ingrowth.
8. Close the incision site using standard surgical technique.
9. Discard any unused portions of the scaffold per your institution's procedures for biohazardous materials.

HOW SUPPLIED

GalaFLEX 3DR™ scaffold is available as a sterile, undyed scaffold of varying widths and lengths.

STERILISATION

GalaFLEX 3DR™ scaffold is sterilised using Ethylene Oxide. GalaFLEX 3DR™ scaffold is supplied sterile and is for SINGLE USE ONLY.

DO NOT clean, re-sterilise or re-use as this may damage or compromise performance of devices and may expose patient to risk of transmitting infectious disease.

STORAGE

Store at room temperature, 15° to 25°C. Avoid prolonged exposure to elevated temperatures.

Printed on: 10 Apr 2024, 11:37:49 am; Printed by: KEARLE.

SYMBOLS KEY

	Catalogue Number
	Batch Code
	Use By – year, month and day
	Do Not Re-use
	Do not use if package is damaged
	Temperature limitation
R_x Only	Rx ONLY Prescription Use Only
	Sterilised by Ethylene Oxide
	Do Not Re-sterilise
	CAUTION - Consult Instructions for Use
	Manufacturer
	European Union Authorised Representative
	CE mark and identification number of Notified Body
	MR Safe

CONTACT INFORMATION

Manufacturer

Tepha, Inc.
99 Hayden Avenue, Suite 360
Lexington, MA 02421 USA

Galatea Surgical, Inc. is a wholly owned subsidiary of Tepha, Inc.

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Printed on: 10 Apr 2024, 11:37:49 am; Printed by: KEARLE.

Released

Always confirm use of the current revision of this document. 400228LX, 0

www.galateasurgicalinternational.com

Released

400228 R0
2023/05

Printed on: 10 Apr 2024, 11:37:49 am; Printed by: KEARLE.

Signature Manifest

Document Number: 400228LX

Revision: 0

Title: GalaFLEX 3DR Scaffold eIFU, CE

All dates and times are in Eastern Standard Time.

DA CR-21370-EC-GLFLX-IFU

DC Change Request

Name/Signature	Title	Date	Meaning/Reason
Brenda Schultz (10172834)	Coord 2, Quality Document Sys		
Documentation Placeholder (DOCHOLDER)	Place Holder		
Julia Kue (JKUE)	Documentation Lead		
Michelle DAmbra (10300496)	Coord 2, Quality Document Sys		
Serenna Boland (10372452)	QA Doc System Mngmt 2	26 May 2023, 10:47:28 AM	Approved

DC Check In 2

Name/Signature	Title	Date	Meaning/Reason
April Fry (10331175)	Global Labeling Manager		
Philip Gambrell (PGAMBRELL)	Principal AQE		
Tara Gatta (TGATTA)	PKg Sr. Global Lbl Spec.	01 Jun 2023, 10:21:06 AM	Complete

Team Review/Approval

Name/Signature	Title	Date	Meaning/Reason
April Fry (10331175)	Global Labeling Manager	02 Jun 2023, 01:01:18 PM	Approved
Alyssa Watson (AWATSON)	Sr. Mgr Quality	02 Jun 2023, 02:17:53 PM	Approved
Adam Cavanagh (ACAVANAGH)	Sr Program Mgr-Global MKT Mang	08 Jun 2023, 04:22:26 PM	Approved
Adam Mercer (10367662)	Staff Regulatory Affairs Spe	28 Jul 2023, 10:19:51 AM	Approved
Tara Gatta (TGATTA)	PKg Sr. Global Lbl Spec.	28 Jul 2023, 10:32:59 AM	Approved

RA Attachments

Name/Signature	Title	Date	Meaning/Reason
Adam Mercer (10367662)	Staff Regulatory Affairs Spe	28 Jul 2023, 11:44:53 AM	Complete

DC Final Review

Name/Signature	Title	Date	Meaning/Reason
Brenda Schultz (10172834)	Coord 2, Quality Document Sys		
Documentation Placeholder (DOCHOLDER)	Place Holder		
Michelle DAmbra (10300496)	Coord 2, Quality Document Sys	28 Jul 2023, 12:40:01 PM	Complete

Release

Name/Signature	Title	Date	Meaning/Reason

Brenda Schultz (10172834)	Coord 2, Quality Document Sys	28 Jul 2023, 12:46:38 PM	Approved
---------------------------	-------------------------------	--------------------------	----------

Notify

Name/Signature	Title	Date	Meaning/Reason
Tony Rennick (10308044)	Pkg Sr Global Lbl Spec	28 Jul 2023, 12:46:38 PM	Email Sent
Melissa Liquori (10203502)	Sr Quality Engineer	28 Jul 2023, 12:46:38 PM	Email Sent
Skander Limem (10203501)	Director of New Product Innova	28 Jul 2023, 12:46:38 PM	Email Sent
Tara Gatta (TGATTA)	PKg Sr. Global Lbl Spec.	28 Jul 2023, 12:46:38 PM	Email Sent

DA CR-21682 admin

DC Change Request

Name/Signature	Title	Date	Meaning/Reason
Brenda Schultz (10172834)	Coord 2, Quality Document Sys		
Documentation Placeholder (DOCHOLDER)	Place Holder		
Julia Kue (JKUE)	Documentation Lead		
Serenna Boland (10372452)	QA Doc System Mngmt 2		
Jeremy Motta (10203510)	Sr. Quality System Eng		
Michelle DAmBra (10300496)	Coord 2, Quality Document Sys	02 Aug 2023, 10:48:19 AM	Approved

DC Check In 2

Name/Signature	Title	Date	Meaning/Reason
Documentation Placeholder (DOCHOLDER)	Place Holder		
Michelle DAmBra (10300496)	Coord 2, Quality Document Sys	02 Aug 2023, 10:53:44 AM	Complete

Originator Approval

Name/Signature	Title	Date	Meaning/Reason
Michelle DAmBra (10300496)	Coord 2, Quality Document Sys	02 Aug 2023, 11:01:31 AM	Approved
Adam Mercer (10367662)	Staff Regulatory Affairs Spe	02 Aug 2023, 05:14:02 PM	Approved

DC Final Release

Name/Signature	Title	Date	Meaning/Reason
Brenda Schultz (10172834)	Coord 2, Quality Document Sys		
Documentation Placeholder (DOCHOLDER)	Place Holder		
Julia Kue (JKUE)	Documentation Lead		
Serenna Boland (10372452)	QA Doc System Mngmt 2		
Jeremy Motta (10203510)	Sr. Quality System Eng		
Michelle DAmBra (10300496)	Coord 2, Quality Document Sys	03 Aug 2023, 07:44:27 AM	Approved

Notify

Name/Signature	Title	Date	Meaning/Reason
Michelle DAmBra (10300496)	Coord 2, Quality Document Sys	03 Aug 2023, 07:44:27 AM	Email Sent