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**ENGLISH**

**DESCRIPTION**

GalaFLEX® scaffold is a single layer, bioresorbable surgical mesh made from poly-4-hydroxybutyrate (P4HB). P4HB is produced from a naturally occurring monomer and is processed into monofilament fibres and knitted into a surgical scaffold. The structure of the scaffold consists of a knit pattern that mitigates the potential for unravelling when the scaffold is cut. The scaffold was fabricated specifically to provide strength throughout the critical wound healing period in soft tissue repair and reinforcement indications.

The design of the scaffold allows for an open pore pattern that encourages the proliferation of healthy tissue at the repair site, after which the product completely bioresorbs. Although the scaffold loses strength with time, its porous construction was designed to allow native tissue ingrowth and gradual transfer of load from the scaffold to the tissue.

**INTENDED USE**

GalaFLEX® scaffold is intended for use, as an adjunct to sutures, for the reinforcement and repair of soft tissue where weakness exists and where the addition of a reinforcing material is needed to obtain the desired surgical result in patients undergoing breast surgery.

The GalaFLEX® scaffold is designed to be used in patients undergoing soft tissue repair and reinforcement in medically necessary breast surgery procedures where the existing soft tissue is deficient to support the surgical repair. Examples of such breast surgery applications include reduction mammoplasty and breast revision surgery to correct a medical condition. GalaFLEX® scaffold may also be used in cosmetic breast procedures.

**CONTRAINDICATIONS**

None known.

**WARNINGS**

1. Device manufacture involves exposure to tetracycline hydrochloride and kanamycin sulfate. The safety and product use for patients with hypersensitivities to these antibiotics is unknown.
2. The safety and effectiveness of GalaFLEX® scaffold in neural tissue and in cardiovascular tissue has not been established.
3. The safety and effectiveness of GalaFLEX® scaffold in paediatric use has not been established.
4. Placement of the scaffold in direct contact with bowel or viscera is not recommended.
5. If an infection develops, treat the infection aggressively. An unresolved infection may require removal of the scaffold.
6. Because GalaFLEX® scaffold is fully bioresorbable, it should not be used in repairs where permanent support from the scaffold is required.
7. GalaFLEX® scaffold is supplied sterile. Inspect the device and packaging prior to use to be sure they are intact and undamaged.

8. GalaFLEX® is for single use only. Do not re-sterilise or re-use any portion of GalaFLEX® scaffold. Never re-use a scaffold, even if it seems to be in perfect condition, in order to prevent risk of cross-contamination or risk of reduced performance. Do not re-sterilise the scaffold. Once the package has been opened, any unused scaffold should be discarded.
9. Unused scaffold must be discarded according to the institution's procedures for handling of biohazardous materials.
10. Any decision to remove the scaffold should take into account potential risks associated with a second surgical procedure. This may include difficulty to remove the scaffold due to the ingrowth of tissue or degradation of scaffold. Scaffold removal should be followed by adequate post-operative management.
11. The GalaFLEX® scaffold has not been studied for use in breast reconstructive surgeries.

#### **PRECAUTIONS**

Only doctors qualified in the appropriate surgical techniques should use this device. Users should be familiar with surgical procedures and techniques, including strength requirements and scaffold size choices. Improper selection, placement, positioning and fixation of GalaFLEX® can cause subsequent undesirable results. Patient should be provided with instructions regarding post-operative care, i.e. lifting, hygiene, activity limitations and any other specific patient/procedure requirements. GalaFLEX® scaffold is bioresorbable and therefore any initial palpability of the scaffold decreases over time.

The safety and effectiveness of GalaFLEX® scaffold in the proximity of existing or excised cancer has not been established.

#### **ACTIONS**

GalaFLEX® scaffold degrades through a process of hydrolysis and hydrolytic enzymatic digestion process.

It has been developed to minimise the variability of resorption rate and strength and provide support throughout the expected period of healing.

Pre-clinical implantation studies indicate that GalaFLEX® scaffold retains approximately 70% of its strength at 12 weeks. Absorption of the scaffold material will be essentially complete within 18-24 months.

#### **ADVERSE REACTIONS**

In pre-clinical testing, GalaFLEX® scaffold elicited a minimal tissue reaction characteristic of foreign body response to a substance. The tissue reaction resolved as the scaffold was resorbed. Possible complications of using GalaFLEX® scaffold include infection, seroma, pain or swelling, scaffold migration, wound dehiscence, haemorrhage, adhesions, haematoma, inflammation, extrusion and recurrence of the soft tissue defect.

#### **IMAGING COMPATIBILITY**

GalaFLEX® scaffold is made from fully resorbable P4HB and is temporary in nature. This material is non-conducting, non-metallic and non-magnetic. Therefore, in accordance with the definition stated in ASTM F-2503-13,

Standard Practice for Marking Devices and Other Items for Safety in the Magnetic Resonance Environment, the device is determined to be "MR Safe - an item that poses no known hazards in all MR environments."

**DIRECTIONS FOR USE**

1. Prior to implanting GalaFLEX® scaffold, complete repair using appropriate suture technique.
2. GalaFLEX® scaffold should be prepared for placement using standard surgical preparation technique.
3. Using aseptic technique, GalaFLEX® scaffold may be cut to the shape or size desired for a specific application.
4. To ensure the scaffold remains flat, minimise pulling and stretching of the scaffold prior to implantation.
5. Transfer the scaffold to the surgical site. Implant the scaffold so that it lays flat against the tissue and the edge of the scaffold extends beyond the margins of the defect. Suture the scaffold into place, avoiding excessive tension and without the expectation of stretch.
6. It is recommended that suture fixation be placed 6 mm to 12 mm apart at a distance approximately 6 mm from the edge of the scaffold using an interrupted suturing technique.
7. If using a resorbable suture to fixate the scaffold, please ensure the strength retention of the suture is sufficient to secure the scaffold during the healing period.
8. Close the incision site using standard surgical technique.
9. Discard any unused portions of the scaffold per your institution's procedures for biohazardous materials.

**HOW SUPPLIED**

GalaFLEX® scaffold is available in single packets as a sterile, undyed scaffold in single sheet sizes of varying widths and lengths.

**STERILISATION**

GalaFLEX® scaffold is sterilised using Ethylene Oxide. GalaFLEX® scaffold is supplied sterile and is for SINGLE USE ONLY.

DO NOT clean, re-sterilise or re-use as this may damage or compromise performance of devices and may expose patient to risk of transmitting infectious disease.

**STORAGE**

Store at room temperature, 15° to 25°C. Avoid prolonged exposure to elevated temperatures.

**SYMBOLS KEY**

	Catalogue Number
	Batch Code
	Use By – year, month and day
	Do Not Re-use
	Do not use if package is damaged
	Temperature limitation
<b>R<sub>x</sub> Only</b>	Rx ONLY Prescription Use Only
	Sterilised by Ethylene Oxide
	Do Not Re-sterilise
	CAUTION - Consult Instructions for Use
	Manufacturer
	European Union Authorised Representative
	CE mark and identification number of Notified Body
	MR Safe

**CONTACT INFORMATION**

Manufacturer

Tepha, Inc.  
99 Hayden Avenue, Suite 360  
Lexington, MA 02421 USA

Galatea Surgical, Inc. is a wholly owned  
subsidiary of Tepha, Inc.

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Printed on: 23 Dec 2024, 12:32:33 pm; Printed by: KEARLE.

Released

Always confirm use of the current revision of this document. 400179LX, 1

[www.galateasurgicalinternational.com](http://www.galateasurgicalinternational.com)

Released

400179 R1  
2024/04

Printed on: 23 Dec 2024, 12:32:33 pm; Printed by: KEARLE.

## Signature Manifest

**Document Number:** 400179LX

**Revision:** 1

**Title:** GalaFLEX Scaffold eIFU, CE

All dates and times are in Eastern Standard Time.

### DA CR-22614 (BT21 Relaunch)-LX1

#### DC Change Request

Name/Signature	Title	Date	Meaning/Reason
Documentation Placeholder (DOCHOLDER)	Place Holder		
Michelle DAmbra (10300496)	Coord 2, Quality Document Sys	04 Apr 2024, 08:45:51 AM	Approved

#### Team Review/Approval

Name/Signature	Title	Date	Meaning/Reason
Melissa Liquori (10203502)	Sr Quality Engineer	10 Apr 2024, 08:33:46 AM	Approved
James McAdam (JMCADAM)	PKG ENGINEER, PACKAGING	10 Apr 2024, 12:40:36 PM	Approved
Katherine Earle (KEARLE)	Mgr, Regulatory Affairs	10 Apr 2024, 02:59:33 PM	Approved
Tara Gatta (TGATTA)	PKg Sr. Global Lbl Spec.	10 Apr 2024, 04:47:20 PM	Approved
Martine Rizzo (10358800)	Associate Product Manager	10 Apr 2024, 08:54:04 PM	Approved
April Fry (10331175)	Global Labeling Manager	11 Apr 2024, 11:35:04 AM	Approved

#### Release

Name/Signature	Title	Date	Meaning/Reason
Serenna Boland (10372452)	QA Doc System Mngmt 2	11 Apr 2024, 01:36:20 PM	Approved

#### Notify

Name/Signature	Title	Date	Meaning/Reason
Adam Mercer (10367662)	Staff Regulatory Affairs Spe	11 Apr 2024, 01:36:21 PM	Email Sent
Melissa Liquori (10203502)	Sr Quality Engineer	11 Apr 2024, 01:36:21 PM	Email Sent
April Fry (10331175)	Global Labeling Manager	11 Apr 2024, 01:36:21 PM	Email Sent